


Bestrahlungsdauer in Abhängigkeit der Schalterstellungen

in (Stunden :) Minuten : Sekunden

Modell CL plus 100:



Dosis (J/cm ²): ►	0,5	1	2	3	4	5	6	8	
									Leistung:
Akup. Kind	00:05	00:05	00:05	00:05	00:05	00:05	00:05	00:05	25%
Akup. Erw.	00:15	00:15	00:15	00:15	00:15	00:15	00:15	00:15	50%
1 cm ²	00:05	00:10	00:20	00:30	00:40	00:50	01:00	01:20	100%
2 cm ²	00:10	00:20	00:40	01:00	01:20	01:40	02:00	02:40	100%
5 cm ²	00:25	00:50	01:40	02:30	03:20	04:10	05:00	06:40	100%
10 cm ²	00:50	01:40	03:20	05:00	06:40	08:20	10:00	13:20	100%
20 cm ²	01:40	03:20	06:40	10:00	13:20	16:40	20:00	26:40	100%
30 cm ²	02:30	05:00	10:00	15:00	20:00	25:00	30:00	40:00	100%
50 cm ²	04:10	08:20	16:40	25:00	33:20	41:40	50:00	01:06:40	100%